OstiumMAX®
OstiumMAXG®
OstiumMAXBlock®

ESP

IMPLANTE DE MATRIZ ÓSEA BOVINA

PRESENTACIONES:

OstiumMAX

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Vial monodosis de 0,5g

Vial monodosis de 1g

OstiumMAX G

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Vial monodosis de 1g

Vial monodosis de 1,5g

Vial monodosis de 2g

OstiumMAX Block

Estructuras óseas

COMPOSICIÓN:

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Bone Structures

COMPOSITION:

OstiumMAX is a mineral from bovine bone. The highly purified mineral osteoconductive structure is made from natural bone through a multi-stage purification process that strictly complies with safety standards. Due to its natural origin, OstiumMAX is comparable to mineralized human bone (natural nanocrystalline apatite), from both the chemical and structural points of view.

PROPERTIES / EFFECTS:

THE PRODUCT IS PROVIDED STERILE. DO NOT RE-STERILIZE. NONTOXIC

The inorganic bone matrix of Ostium Max has a macroporous and microporous structure similar to the human cancellous bone. The natural composition and the large volume of interconnected pores encourage the formation and growth of new bone at the implantation site. Over time, Ostium Max undergoes partial remodeling by the action of osteoclasts and osteoblasts (physiological remodeling). Due to its properties, Ostium Max is a viable alternative to autogenous bone in appropriate defects.

INDICATION:

OstiumMAX is indicated as a bone substitute and/or bone filling with osteoconductive and osteoinductive action.

- Augmentation/reconstruction of the alveolar ridge;
- Post-extraction socket filling;
- Implantology: Implant site preparation, filling of bone dehiscences and maxillary sinus floor elevation;
- Periodontology: filling of bone defects, membrane holder during guided tissue regeneration (GTR).

Ostium Max should only be used by authorized professionals trained in this field as dentists and surgeons.

INSTRUCTIONS FOR USE:

To use Ostium Max, the general guidelines for sterile handling and patient medication should be followed.

- Complete elimination of granulation tissue after defect exposure.
- Mixture of Ostium Max granulated with the patient's blood or with saline solution prior to implantation.

APPLICATION:

OstiumMax is applied to the defect using sterile instruments (forceps, spatula, spoon or syringe).

It is advisable to cover Ostium Max with a membrane barrier (OstiumMAX Cover) to act as mechanical support and promote adhesion. When closing the wound, the flap of soft tissue should completely cover the Ostium Max implant and should be fixed with sutures.

If the primary closure of the wound is not complete, the flap should be mobilized (by an incision through the periosteum or covering the wound with a membrane).

SPECIAL INSTRUCTIONS FOR USE IN PERIODONTOLOGY:

A hygiene phase is recommended before surgery, which will include the provision of appropriate instructions to the patient, since a basic requirement for the success of the treatment is the control of bacterial infections as a meticulous oral hygiene. A postoperative maintenance phase can ensure a long-term successful treatment. In addition to plaque control, periodontal defect fill with Ostium-Max requires a successful local treatment of the periodontal lesion

(root planing, debridement) before implantation. To achieve optimal tissue regeneration, the defect should be covered with a membrane.

LIMITATIONS OF USE/PRECAUTIONS:

Contraindications:

OstiumMAX must not be used in infected areas.

Precautions:

Material for single use.

Material for use in its entirety on a single receiver.

Its reuse is prohibited.

Do not use if package is not complete.

Sterile. Its reesterilization is prohibited.

Read the instructions before using the product.

OstiumMAX must be used with special caution in patients with:

- acute or chronic infection at the surgical site. (E.g., osteomyelitis)
- uncontrolled metabolic diseases such as diabetes, osteomalacia, severe renal or liver disease, thyroid disorder
- auto immunediseases
- radiotherapy
- prolonged corticosteroid therapy
- high consumption of tobacco or alcohol
- malnutrition
- suspected foreign body sensitivity or allergies. (Studies must be done to determine what material to implant)
- psychological condition that prevents treatment

It should only be implemented in vital bone tissue in direct contact with the host bone (if necessary, with microfracture of the bone surface) to ensure bone regeneration. For larger defects, the addition of autogenous cancellous bone or bone marrow can improve the regeneration process. Experience shows that mobilization by mechanical loads (compression loads) or the insertion of implants (two-stage procedure) in the enlarged area must be avoided for at least 6 months after the insertion of OstiumMAX.

Pregnancy / Breastfeeding:

There is no data available on the use of the product during pregnancy or breastfeeding. For safety reasons: pregnant or nursing women must not be treated with OstiumMAX. The safety and efficacy of Ostium-

MAX have not been researched in children before skeletal maturity.

SIDE EFFECTS:

NO ACCIDENTS OR ADVERSE EFFECTS SHOULD OCCUR IF THE INSTRUCTIONS AND PRECAUTIONS DESCRIBED ABOVE ARE FOLLOWED.

Some of the undesired effects may include:

- infection
- persistent pain
- failure of osseointegration
- nerve damage

STORE IN A DRY PLACE AT ROOM TEMPERATURE

Do not use after the expiration date on the container.

Sterilized by: Gamma radiation.

Technical Director: Pharmacist, Lattini Natacha Verónica, License No. 3838

Condition of sale: Exclusive sale to professionals and health institutions.

Authorized by ANMAT PM - 2158-1

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